


PROFEEL[®]

DHD[™] EXTRA PROTECTION

Powder Free Latex Surgical Gloves



Damp Hand Donning - DHD[™] • Natural Rubber Latex • Powder Free

www.wrpworld.com

CE
2797
BSI (2797)
1066 EP, NL

EN ISO 374-1:2016
Type B
KPT

EN ISO 374-5:2016
VIRUS

PROFEEL[®]

DHD[™] EXTRA PROTECTION

Powder Free Latex Surgical Gloves

- Innovatively designed in understanding need for exceptional protection along with superior tactile sensitivity and comfort.
- Ideal for use during intense orthopedic surgery and high risk situations which require gloves with sufficient thickness and resilience to ensure barrier integrity is not compromised.
- An optimum strength-to-thickness ratio provides appropriate tactile sensitivity to promote dexterity and enhanced precision.
- Damp Hand Donning (DHD)- The inner surface of the glove is coated with a proprietary polymer to enhance damp hand donning.
- Beaded cuff prevents tearing while donning and helps prevent roll back

PRODUCT INFORMATION

Article No.	:	WRP SG NR PF 125
Style No.	:	WRP SG NR PF 125 - EP DHD
Material	:	Natural Rubber Latex
Coating	:	Damp Hand Donning (DHD)
Weight per pair, size 7.5, g	:	23.4 ± 0.6
Sizes	:	5.5 - 9.0
Palm Thickness, mm (min)	:	0.25
Length (mm)	:	290 ± 10
Colours	:	Brown

PRODUCT CONFORMANCE:

Medical Device: in compliance with European Medical Device Directive 93/42/EEC (CE Class IIa)

QUALITY ASSURANCE:

Manufacturing process is in compliance with US FDA Quality System Regulation (QSR), ISO 9001 Quality Management System and BS EN ISO 13485 Quality Systems



 Your Partner In Protection[™]

EC REP

REMESCO Handelsges.m.b.H.

Grinzing Allee 5/25, 1190 Vienna,
AUSTRIA

Office : +43-1-328-50-88

Fax : +43-1-328-29-49

Email : office@remesco.com

WRP Asia Pacific Sdn Bhd 147817V

Lot 1, Jalan 3, Kawasan Perusahaan
Bandar Baru Salak Tinggi,
43900 Sepang, Selangor Darul Ehsan,
MALAYSIA

Telephone: +60-3-8706 1486

Fax : +60-3-8706 1557

Email : customer_wrp@wrpworld.com
sales@wrpworld.com

Website : www.wrpworld.com