

# **PROFEEL<sup>®</sup>**

**DHD<sup>™</sup> Sensitive Powder Free**  
**Aloe Coated Latex Surgical Gloves**

## **Powder Free Latex Sterile Surgical Gloves**

PROFEEL Powder Free  
PROFEEL DHD<sup>™</sup> Platinum  
PROFEEL DHD<sup>™</sup> Micro  
PROFEEL DHD<sup>™</sup> Extra Protection  
PROFEEL DHD<sup>™</sup> Sensitive



Damp Hand Donning – DHD<sup>™</sup> • Natural Rubber Latex • Powder Free



Your Partner In Protection<sup>™</sup>



## Exceptional Dexterity

**PROFEEL<sup>®</sup> DHD™ Sensitive Powder Free Latex Surgical Gloves** are designed with a thin latex film that improves tactile sensitivity and comfort during surgery. Its optimal thickness and unique surface finish maximises sensitivity even when wearing two pairs of gloves.

The combination of sensitivity and Aloe Vera coating, make **PROFEEL<sup>®</sup> DHD™ Sensitive Powder Free Latex Surgical Gloves** an ideal underglove.

The special Aloe Vera gel that coats the inside of **PROFEEL<sup>®</sup> DHD™ Sensitive Powder Free Latex Surgical Gloves** naturally moisturises and hydrates your skin, helping to keep it soft and healthy.

- The inner surface is coated with polymer to enhance damp hand donning
- Thinner than normal surgical gloves for added sensitivity
- Aloe Vera coating on inner surface helps moisturize and hydrate the hand
- Anatomically shaped to reduce hand fatigue during prolonged surgical work
- Beaded cuff prevents tearing while donning

## PRODUCT INFORMATION

Article No.	:	WRP SG NR PF 125
Style No.	:	WRP SG NR PF 125 - SN
Material	:	Natural Rubber Latex
Coating	:	Damp Hand Donning (DHD)
Weight per pair, size 7.5, g	:	20.4 ± 0.6
Sizes	:	5.5 - 9.0
Palm Thickness, mm (min)	:	0.15
Length (mm)	:	290 ± 10
Colours	:	Green

## PRODUCT CONFORMANCE:

Medical Device: in compliance with  
European Medical Device Directive  
93/42/EEC (CE Class IIa)

## QUALITY ASSURANCE:

Manufacturing process is in compliance  
with US FDA Quality System Regulation  
(QSR), ISO 9001 Quality Management System  
and BS EN ISO 13485 Quality Systems



### WRP Asia Pacific Sdn Bhd 1 4 7 8 1 7 V

Lot 1, Jalan 3, Kawasan Perusahaan  
Bandar Baru Salak Tinggi,  
43900 Sepang, Selangor Darul Ehsan,  
MALAYSIA

Office +60-3-8706 1486  
Facsimile +60-3-8706 1557  
Email [customer\\_wrp@wrpworld.com](mailto:customer_wrp@wrpworld.com)



### REMESCO Handelges.m.b.H,

Grinzinger Allee 5/25, 1190 Vienna,  
AUSTRIA

Office +43-1-328-50-88  
Facsimile +43-1-328 29-49  
Email [office@remesco.com](mailto:office@remesco.com)

Distributed by: