

PROFEEL[®]

DHD[™] Micro Powder Free
Latex Surgical Gloves



Damp Hand Donning - DHD[™] • Natural Rubber Latex • Powder Free

www.wrpworld.com

CE
2797
BSI (2797)
1066 EP, NL

EN ISO 374-1:2016
Type B

KPT

EN ISO 374-5:2016

VIRUS

PROFEEL[®]

DHD[™] Micro Powder Free

Latex Surgical Gloves

- Specially designed for specialists who demand extra sensitivity in the fingertips for micro and ophthalmic surgeries.
- With WRP's new PROFEEL[®] DHD[™] Polyisoprene Surgical Gloves, you can enjoy the fit, feel and comfort of natural rubber latex in a 100% latex free product.
- Powder free feature eliminates powder-induced irritation and powder stains
- Environmentally friendly and biodegradable
- A proprietary manufacturing process that ensures the exterior of our gloves consistently match the smooth grip of the human hand
- Damp Hand Donning (DHD)- The inner surface of the glove is coated with a proprietary polymer to enhance damp hand donning
- Beaded cuff prevents tearing while donning and helps prevent roll back

PRODUCT INFORMATION

Article No.	: WRP SG NR PF 125
Style No.	: WRP SG NR PF 125 - MICRO DHD
Material	: Natural Rubber Latex
Coating	: Damp Hand Donning (DHD)
Weight per pair, size 7.5, g	: 23.4 ± 0.6
Sizes	: 5.5 - 9.0
Palm Thickness, mm (min)	: 0.15
Length (mm)	: 290 ± 10
Colours	: Brown

Glove Size	Palm Width (mm)
5 ½	72 ± 4
6	77 ± 5
6 ½	83 ± 5
7	89 ± 5
7 ½	95 ± 5
8	102 ± 6
8 ½	108 ± 6
9	114 ± 6

PRODUCT CONFORMANCE:

Medical Device: In compliance with European Medical Device Directive 93/42/EEC (Class IIa), EN 455: Part 1/2/3/4, EU Regulation 2016/425, EN 420:2003 +A1 2009, EN ISO 374-1:2016, EN ISO 374-2:2014, EN ISO 374-4:2013, EN 374-5:2016, EN 16523-1:2015, ASTM D3577, FDA 510 K

QUALITY ASSURANCE:

Manufacturing process is in compliance with US FDA Quality System Regulation (QSR), ISO 9001 Quality Management System and BS EN ISO 13485 Quality Systems



 Your Partner In Protection[™]

EC REP

REMESCO Handelsges.m.b.H.

Grinzinger Allee 5/25, 1190 Vienna,
AUSTRIA

Office : +43-1-328-50-88

Fax : +43-1-328-29-49

Email : office@remesco.com

WRP Asia Pacific Sdn Bhd 147817V

Lot 1, Jalan 3, Kawasan Perusahaan
Bandar Baru Salak Tinggi,
43900 Sepang, Selangor Darul Ehsan,
MALAYSIA

Telephone : +60-3-8706 1486

Fax : +60-3-8706 1557

Email : customer_wrp@wrpworld.com
sales@wrpworld.com

Website : www.wrpworld.com