

# PROFEEL®

## DOUBLE GLOVING SYSTEM

FOR PUNCTURE DETECTION  
DHD LATEX SURGICAL UNDERGLOVE AND OVERGLOVE

Powder Free Sterile Surgical Gloves

PROFEEL Powder Free

PROFEEL DHD Platinum

PROFEEL DHD Micro

PROFEEL DHD Extra Protection

PROFEEL DHD Sensitive

PROFEEL DHD Double Gloving System




Your Partner In Protection™



# PROFEEL®

## DOUBLE GLOVING SYSTEM

FOR PUNCTURE DETECTION  
DHD LATEX SURGICAL UNDERGLOVE AND OVERGLOVE



### Exceptional Protection

- Two pairs of gloves per pack: overgloves and undergloves
- Wearing two pairs of surgical gloves provide an additional barrier and further reduce the risk of contamination.
- Reduce the number of perforation to innermost pair of surgical gloves.
- Natural Rubber Latex's intrinsic softness provides an exceptionally comfortable and flexible glove.
- Damp Hand Donning (DHD™) – Inner surface of the glove is coated with a proprietary polymer to enhance damp hand donning.
- A non slick skin-grip surface finish provides excellent instrument grip, under both wet and dry conditions.
- A superior tactile sensitivity provides optimum 'feel'.
- Specially designed high tack band near the cuff of the glove to securely hold onto the wearer's gown to minimize roll down.
- A reinforced beaded cuff feature ensures easy donning while preventing tearing.
- Tested for water leak before packing to AQL 0.65.
- The cuff is printed with size and logo for easy identification.

### Dimensions (Overglove & Underglove)

Glove Size	Reorder No.	Palm Width (mm)	Length (mm)
5½	P3355-61	72±4	Min 278
6	P3360-61	77±5	Min 280
6½	P3365-61	83±5	Min 280
7	P3370-61	89±5	Min 283
7½	P3375-61	95±5	Min 287
8	P3380-61	102±6	Min 288
8½	P3385-61	108±6	Min 290

### Typical Thickness

Location of Thickness Measurements	Overglove	Underglove
Finger	0.21 ± 0.02	0.18 ± 0.03
Palm	0.17 ± 0.02	min. 0.10
Cuff	0.16 ± 0.02	min. 0.10

### Typical Physical Properties

Parameters	Overglove		Underglove	
	Before Aging	After Aging	Before Aging	After Aging
Tensile Strength (MPa)	Min. 24	Min. 18	Min. 24	Min. 18
Elongation (%)	Min. 750	Min. 560	Min. 560	Min. 560
Force at Break (N)	Min. 12	Min. 9	Min. 9	Min. 9

### Pre-shipment Quality Inspection

Criteria	Insp. Level	AQL
Dimension	N=13 (EN455-2); S-2	4.0
Physical Properties	N=13 (EN455-2); S-2	4.0
1000ml Water Leak Test	G-I	0.65
Visual Inspection for Major Defects	G-I	1.5
Visual Inspection for Minor Defects	G-I	1.0

### PRODUCT INFORMATION

#### Packaging Information

: One pouch contains one pair of Overgloves as per size printed on the packaging and one pair of Undergloves half a size larger

#### Overglove

#### Type

: Powder free, polymer coated, damp-hand-donnable (DHD™) and Gamma-sterilized surgical gloves

#### Primary Material

: Natural rubber latex

#### Color

: Off white

#### Design and Features

: Hand Specific, curved fingers, bisque finish (SH3) and beaded cuff

#### Surface Treatment

##### Underglove

: Polymer coating on inner surface to ease donning

#### Type

: Micro thin, powder free, polymer coated, damp-hand-donnable (DHD) and Gamma-sterilized surgical gloves

#### Primary Material

: Natural rubber latex

#### Color

: Green (PMS 3275U)

#### Design and Features

: Hand Specific, curved fingers, smooth finish for easy double gloving and beaded cuff

#### Surface Treatment

: Polymer coated on inner surface to ease donning. Outer surface is smooth



**WRP Asia Pacific Sdn Bhd** 147817V

Lot 1, Jalan 3, Kawasan Perusahaan Bandar Baru Salak Tinggi  
43900 Sepang, Selangor Darul Ehsan  
MALAYSIA

Office +60-3-8706 1486  
Facsimile +60-3-8706 1557  
Email customer\_wrp@wrpworld.com



**REMESCO Handelges.m.b.H.**

Grinzinger Allee 5/25, 1190 Vienna,  
AUSTRIA

Office +43-1-328-50-88  
Facsimile +43-1-328 29-49  
Email office@remesco.com

Distributed by: