

trackman®

CONNECTED

THE KIT TO CONNECT YOUR BENCH

Pipette Smarter. Report Faster. Confidently.

I Boost Your Experiment Efficiency

TRACKMAN® Connected guides you through your microplate method to prevent errors. When paired with PIPETMAN® M Connected*, the app displays where, when, and how much to pipette with real time interactions. For more efficiency, verify your pipettes and check your pipette status with the PipetteScope® app.

I Accelerate Collaborative Reporting

- Easily create and share your protocols.
- Export your protocol steps to save locally or for sharing.
- Sync your reports with SciNote Electronic Lab Notebook.

I Improve Traceability on the Bench

Integrate a tablet on the bench in a safe and convenient way. Easily capture notes, pictures, and monitor your protocols. Automatically track protocol run conditions from PIPETMAN M Connected and the connected environment sensor.



*sold separately.

 **GILSON®**



Inside Your Kit



Tablet with PipettePilot and PipetteScope apps

Interact with all your apps and the web at the bench with the most recent SAMSUNG tab.

The tablet comes preloaded with PipettePilot to automatically create and track microplate pipetting tasks and monitor other device activities.

Also loaded on the tablet is PipetteScope, an app that keeps tabs on your pipette's last calibration, user's guides, and other useful pipetting metrics.

Tablet Holder

Microplate Adapter



Strip Adapter

Microplate Holder

Benchtop Pipetting Accessories

Get all the accessories to bring your tablet on the bench and to pipette with it: microplate holder, tempered glass screen protector, 5V power adapter, microplate and strip adapters.

Environment Sensor

Monitor temperature, pressure, and humidity during your experiment with your environment sensor.



Environment Sensor

| TRACKMAN® Connected | Ordering Information |
|----------------------------|----------------------|
| TRACKMAN Connected EU Plug | FB1020 |
| TRACKMAN Connected UK Plug | FB1021 |
| TRACKMAN Connected US Plug | FB1022 |

Make Pipetting Easier for Your Molecular Biology Applications

We've worked closely with QIAGEN to program their gold-standard protocols into TRACKMAN Connected to reduce errors and bring greater traceability to your molecular biology applications.

Depending on your protocol needs, we have three TRACKMAN Connected Starter Kits available: Single Channel, Multichannel, and Pro. Each starter kit includes PIPETMAN M Connected pipettes, TRACKMAN Connected tablet and accessories, and tips.



Learn more at gilson.com/gilson-connect/qiagen

| TRACKMAN® Connected Starter Kit | Ordering Information | |
|---|----------------------|---------|
| TRACKMAN Connected Starter Kit Single Channel | EU | FB10200 |
| | UK | FB10201 |
| | US | FB10202 |
| TRACKMAN Connected Starter Kit Multichannel | EU | FB10210 |
| | UK | FB10211 |
| | US | FB10212 |
| TRACKMAN Connected Starter Kit Pro | EU | FB10220 |
| | UK | FB10221 |
| | US | FB10222 |