

SANDEL[®] Respiratory face masks

M E D I C A L F A C E M A S K S

Protecting healthcare staff and patients from airborne diseases

Protecting healthcare staff and patients from airborne diseases
Respirators should be worn by frontline staff when carrying out a potentially infectious aerosol-generating procedure, when a patient is known or suspected to have an infection spread via the aerosol route or, when caring for patients known or suspected to be infected with a newly identified respiratory virus.

Respirator features include:

- Duckbill style with headbands
- 3 levels of protection : FFP1, FFP2, FFP3
- Resistant to fluid projection (≥ 21.3 kPa or 160mmHg, as per ISO 22609 requirements)
- 2 sizes of FFP2 mask to offer the best fit to wearer's face

KEY FEATURES & BENEFITS

- Duckbill style with headbands
- 3 levels of protection : FFP1, FFP2, FFP3
- Resistant to fluid projection (≥ 21.3 kPa or 160mmHg, as per ISO 22609 requirements)
- 2 sizes of FFP2 mask to offer the best fit to wearer's face

RECOMMENDED FOR

- Enhancing infection control
- Preventing the risk of cross-contamination



SANDEL[®] Anticipate the Unexpected

Innovative SANDEL[®] Safety Solutions are developed through collaboration with clinical professionals to help prevent injuries, infections, and errors in high-risk healthcare environments. SANDEL products provide the confidence that healthcare professionals need to perform at their best and protect patients from preventable harm.

SANDEL[®] Respiratory face masks

M E D I C A L F A C E M A S K S

PRODUCT DESCRIPTION

Color	Orange, Green, White
Shape	Duckbill
String Type	Elastic bands

MANUFACTURING AND SAFETY STANDARDS

Lot Traceability	Yes
Reusable	No, Disposable Product, Single Use Only
Resterilization	Not Applicable
Shape	Duckbill
String Type	Elastic bands
Standards	ISO 9001, EN 149:2001 + A1:2009, ISO 22609:2004

PACKAGING AND STORAGE

Packaging Overview	50 pcs / box 10 boxes / case
Case Quantity	500, 500, 500, 500
Shelf Life	5 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition

Contact your Ansell representative for ordering or more information.

North America

Tel: (866) 764-3327
Email: sandel-cs@ansell.com
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean

Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa

Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific

Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand

Tel: +61 3 9270 7270
Email: protection@ap.ansell.com