



PRODUCT DESCRIPTION

Material	Nitrile
Colour	Pink
Shape	Ambidextrous
Cuff	Beaded
External Surface	Finger textured, Chlorinated
Internal Surface	Chlorinated

SIZE / PRODUCT CODE

small	4511	medium	4512
large	4513	x-large	4514

MICRO-TOUCH® NITRAFREE™

A COMBINATION OF ALLERGY PREVENTION AND COMFORT: EVERY HEALTHCARE PROFESSIONAL DESERVES OPTIMAL PROTECTION

- Accelerator free nitrile gloves featuring a soft formulation for increased comfort levels
- Pink coloured nitrile glove for easy identification as a synthetic glove
- High resistance to punctures and tears
- Excellent resistance to a broad range of commonly used chemicals including chemotherapy drugs*
- Textured fingers for excellent wet and dry grip
- Accelerator Free to minimize the risks associated with Type IV allergic reactions
- Beaded cuff
- Both powder and latex-free
- Ambidextrous
- Layer pack for easier dispensing and reduction of waste - save on glove wastage.

Micro-Touch® NitraFree™ pink gloves are both latex and accelerator-free making them a unique examination glove and the product of choice for those suffering from latex or chemical accelerators sensitivity, the latter of which can be confused with latex sensitivity in the absence of testing. Micro-Touch® NitraFree™ also features a softer nitrile formulation for increased comfort levels while maintaining nitrile's characteristic high resistance to tear, puncture and chemicals, including chemotherapy drugs*.

Nitrile Powder-Free Non-Sterile Examination Glove

MD Class I
If relevant, read carefully the instructions for use flyer accompanying the product. Single use only.

Ref-XXXX-20200225

Ansell is a world leader in providing superior health and safety protection solutions that enhance human well being. With operations in North America, Latin America, EMEA and Asia, Ansell employs more than 13,000 people worldwide and holds leading positions in the industrial and medical gloves market, as well as in the sexual health and wellbeing category worldwide. Ansell operates in four main business segments: Medical Solutions, Industrial Solutions, Specialty Markets and Sexual Wellness. Information on Ansell and its products can be found at www.ansell.eu.

Ansell Healthcare Europe N.V. (European Head Office)

Riverside Business Park, Block J
Boulevard International 55, 1070 Brussels, Belgium
Tel. +32 (0) 2 528 74 00 • Fax +32 (0) 2 528 74 01 • Fax Customer Service +32 (0) 2 528 74 03
<http://www.ansell.eu> • E-mail info@ansell.eu