

MICRO-TOUCH[®] NITRAFREE[™]

Performance Series

Nitrile, powder-free pink examination glove designed to minimize Type IV chemical allergies

This non-sterile glove is made from a durable and soft nitrile formulation. Provides excellent protection in all risk settings. Textured fingertips of the glove provide extra protection for wet or dry grip and against slippage.

KEY FEATURES AND BENEFITS:

- Protects from latex Type I allergy
- Manufactured without sulfur-based chemical accelerators, helping to protect users from Type IV (chemical) allergies
- Tested for use with chemotherapy drugs*

* See product packaging or contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Tested for use with chemotherapy drugs
- Sustained, high-risk examination procedures
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] NITRAFREE[™]

Performance Series

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Pink	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.14 / 5.5	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS, S, M, L, XL	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested

PACKAGING AND STORAGE

Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition

STANDARDS AND CERTIFICATIONS

ASTM D6319, EN 420:2003 + A1:2009, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 13485, ISO 10993-10:2002, ISO 9001

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com