

Sterile Singles

Synthetic, sterile single packed, powder-free examination glove ideal for brief, sterile applications

Synthetic, sterile single packed, powder-free examination glove ideal for brief, sterile applications

- Packaged as a single sterile glove per pouch
- Packaged on backing paper for easy donning
- Cost effective synthetic alternative for sterile use
- Protects from latex Type I allergy



Industries

- Healthcare

Recommended For

- Sterile procedures outside of the O.R.
- Brief, low-risk examination procedures
- Protection from Type I latex allergy in HCW's or patients



Sterile Singles

Synthetic, sterile single packed,
powder-free examination glove ideal
for brief, sterile applications

TECHNICAL DATA SHEET

PRODUCT INFORMATION

	Sterile Singles
Material	EVA (Ethylene Vinyl Acetate)
Color	Clear
Glove Design	, Smooth
Cuff	Straight
Manufacturing/QMS Audit Standards	ISO 13485:2003
Regulatory/Standards Compliance	Category III, EN 455 Part 1, EN 455 Part 3, EN 420:2003 + A1:2009, ISO 10993-10:2002, ISO 13485
Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition
Country of Origin	Taiwan
User Needs Segment	Essential Series
Available sizes	S, M, L

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	278 / 10.9		ASTM D3767,EN420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 455-1
Palm Thickness (mm/inches)	0.03 / 1.2		ASTM D3767,EN420
Finger Thickness (mm/inches)	0.03 / 1.2		ASTM D3767,EN420
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	N/A	N/A	ASTM D412 & D573
Elongation at Break (%)	≥600	628	ASTM D412
Force at break (N)	3	2	EN 455-2

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
Prima 6, Prima Avenue
Block 3512, Jalan Teknokrat 6
63000 Cyberjaya
Selangor, Malaysia
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
111 Wood Avenue South,
Suite 210
Iselin, NJ 08830, USA
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
Blvd. Bernardo Quintana No. 7001-C,
Q7001 Torre II,
Suites 1304, 1305 y 1306,
Col. Centro Sur, c.p. 76079
Queretaro, Qro. Mexico
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
Level 3,678 Victoria Street,
Richmond, Vic, 3121
Australia
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ahcenn PYC
123610 Poccnr, Mockba
Kpachonpechehckar
Ha6. 12, n.3, op. 1304A
Ten. +7 495 258 13 16

Technology



Performance Standards and Regulatory Compliance



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2019 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

