

MICRO-TOUCH[®] Micro-Thin[™] White Nitrile

E s s e n t i a l S e r i e s

Provides excellent dexterity, along with exceptional comfort, fit and feel

MICRO-TOUCH[®] Micro-Thin[™] White Nitrile is a non-sterile glove made with a thinner formulation to provide excellent tactile sensitivity and dexterity along with exceptional comfort, fit and feel. High quality manufacturing ensures this glove will hold up for a variety of clinical procedures.

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from Type I latex allergy in HCW's or patients
- Brief, low-risk examination procedures
- Environmentally friendly facilities



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Micro-Thin[™] White Nitrile

Essential Series

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	White	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.06 / 2.4
Finger Thickness (mm/mil)	0.08 / 3.1	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Silicone Free	Not Tested	Sterile	No
Antistatic	Not Tested	Double Gloving Recommendation	Yes, as underglove
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	Essential Series

ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	313015090	313015080	313015070	313015100	313015060

PACKAGING AND STORAGE

Packaging	200 gloves per dispenser; 10 dispensers per case; 2000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Arrêté du 09/04/04, ASTM D6319, EN 16523-1, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 13485, ISO 9001

Sensitisers

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com