


# 0.5 – 5 ml Pipette tips

## 702595 | 702600

PRODUCT  
INFORMATION




Dimensions in mm

<b>Volume</b>	0.5 – 5 ml
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	bulk
<b>Quality</b>	Produced under cleanroom conditions.

# 0.5 – 5 ml Pipette tips

## 702605

PRODUCT  
INFORMATION



Dimensions in mm

<b>Volume</b>	0.5 – 5 ml
<b>Material</b>	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
<b>Color</b>	transparent
<b>Packaging</b>	racked
<b>Quality</b>	Produced under cleanroom conditions.