


1 – 10 ml Pipette tips

702603 | 702604

PRODUCT
INFORMATION




Dimensions in mm

Volume	1 – 10 ml
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	bulk
Quality	Produced under cleanroom conditions

1 – 10 ml Pipette tips

702606

PRODUCT
INFORMATION



Dimensions in mm

Volume	1 – 10 ml
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Produced under cleanroom conditions