



**Product Specification**  
**Medicine cup, 30 ml, (LxØ): 40 x 37 mm,**  
**graduated, PP, yellow**

**Page 1**



### Product description

Order number	75.1337.504
Product description	Medicine cup, max. working volume: 30 ml, (LxØ): 40 x 37 mm, Ø opening: 37 mm, yellow, graduated, material: PP, 75 piece(s)/bag

### Product characteristics

Graduation	yes
Volume of work	30 ml

### Size

Sample volume	30 ml
Diameter	37 mm
Diameter of opening	37 mm
Height of product	40 mm
Length of product	40 mm
Length excluding cap	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-02-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Medicine cup, 30 ml, (LxØ): 40 x 37 mm,**  
**graduated, PP, yellow**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product yellow

**Purity & certification**

Product category Medical device  
CE certified CE supplier  
Batched yes

**Packaging**

Minimum order qty. 3750  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 3750  
Piece (s) / outer case 3750  
Piece(s) / pallet 180000  
Depth of case 190 mm  
Width of case 390 mm  
Height of case 390 mm  
Case volume 0.0289 cbm  
Weight of product 0.0015 kg  
Weight of case 6 kg  
EAN of inner box 4038917127745  
EAN case 4042836002052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-02-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**