



Product description

Order number	76.1340.560
Product description	Lid, for 30-ml medicine cup, material: PE, white, collection and storage of urine, 75 piece(s)/bag

Product characteristics

Application	collection and storage of urine
Example of use	cap for urine collection tube
Suitable for	30-ml medicine cup

Material & colours

Product material	Polyethylene (PE)
Colour of product	white
Colour of cap	white

Purity & certification

Product category	Medical device
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	90000
Depth of case	390 mm
Width of case	265 mm
Height of case	590 mm
Case volume	0.061 cbm
Weight of product	0.0007 kg
Weight of case	3.64 kg
EAN of inner box	4038917127806
EAN case	4042836007156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

