



Sol-Care™

Fast Release Blood Collection Tube Holder

- ▶ The **Sol-Care™ Fast Release Blood Collection Tube Holder** is a non-sterile reusable plastic tube holder with safety mechanism. It is used in conjunction with multi-sample needle in routine blood sampling procedures. The used needle then can be released simply from the holder by pushing the release mechanism to protect users from needlestick injuries.

Features:

- A non-sterile, reusable device
- Simple one-handed needle release technique
- Designed to attach and hold a multi-sample needle during blood collection and release the used needle directly into sharps container
- Compatible with all sizes of Sol-M™ Multi-sample Needles
- Not made with natural rubber latex, PVC or DEHP





Sol - Care™

Fast Release Blood Collection Tube Holder

Sol-Care™ Fast Release Blood Collection Tube Holder

Ref. number	Product description	Bag	Case
BCA028	Fast Release Blood Collection Tube Holder	20	1000

Intended use: The Sol-Care™ Fast release blood collection tube holder is a non-sterile reusable device designed to attach and hold multi-sample needle or safety multi-sample needle during blood collection. The holder has a quick release feature which allows the contaminated needle to be disengaged from the holder with a simple one-handed push-button technique.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- The Sol-Care™ Fast release reusable holder must be cleaned after each procedure and between each patient to prevent the risk of cross contamination.
- Sol-Care™ Fast release Holder can be re-activated - maximum of 100 times.
- Do not activate the device during blood collection procedure.
- There is no need to clean the device when re-used for 1 patient during 1 procedure (changing needle during 1 blood collection procedure)
- Device can be re-processed till max 15 cleaning cycles. The specific number of reuse cycles was validated.
- Do not soak in bleach overnight.
- Never clean with the needle attached to the holder.
- Do not autoclave the device, do not use in heat sterilizer.
- If the holder is visibly soiled with blood or there is failure of product integrity, or failure to activate the release button, discard the device immediately into sharps container according to facility protocol.
- Check a functionality of the safety needle release button prior to the connection with new needle.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMEA-BC-670-0224

For more information, please contact:



Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B
00-843 Warsaw, Poland



info.eu@sol-m.com



Sol M

Building
a Healthier
Tomorrow