






Multichannel Models

VistaLab Catalog #		Color	Volume Range	Dispensing Increments	Accuracy	Precision
8 Channel	12 Channel					
1160-0020 1160-0020L*	1060-0020 1060-0020L*		0.5 – 20µL	0.05µL	±1.0% at 20µL ±1.5% at 10µL ±3.0% at 2µL	0.40% at 20µL 0.80% at 10µL 1.80% at 2µL
1160-0125 1160-0125L*	1060-0125 1060-0125L*		2 – 125µL	0.5µL	±0.8% at 125µL ±1.0% at 62.5µL ±3.0% at 12.5µL	0.25% at 125µL 0.30% at 62.5µL 1.25% at 12.5µL
1160-0250 1160-0250L*	1060-0250 1060-0250L*		5 – 250µL	1µL	±0.8% at 250µL ±1.0% at 125µL ±2.5% at 25µL	0.25% at 250µL 0.40% at 125µL 1.00% at 25µL
	1060-0850 1060-0850L*		25 – 850µL	5µL	±1.0% at 850µL ±1.25% at 425µL ±3.5% at 85µL	0.30% at 850µL 0.40% at 425µL 1.00% at 85µL
1160-1250 1160-1250L*			25 – 1250µL	5µL	±0.8% at 1250µL ±1.0% at 625µL ±3.5% at 125µL	0.20% at 1250µL 0.30% at 625µL 0.80% at 125µL

