



MICRONIC TUBE STORAGE RACK RANGE

Comorack-96

Holds 96 Individual 1.40ml Round Bottom Tubes!



Ordering information

| Catalog number | Description | Unit |
|----------------|---|----------|
| MP22504 | Comorack-96 (including cover) (polycarbonate) | 10 racks |

Superior features

- Based on the Society for Laboratory Automations and Screening (ANSI/SLAS) standards for storage racks
- Standard available with Non-coded 1.40ml round bottom tubes
- Secure and safe locking system
- Autoclavable up to 100 times
- Easily opened with one hand
- Open rack bottom
- Class 7 cleanroom production
- RNase/DNase and Pyrogen free



30 years of experience
in traceable sample
storage

Class 7 clean room
injection molding
and assembly

Production in
the Netherlands and
the United States